










SURFACES

SPANAmerica

Innovative Solutions.

Surfaces Selection Guide

To be used only as a guide. Final product selection should be based on clinical judgement.

	Comfort and Prevention • Low to moderate risk Braden 13-26	Prevention and Low Level Treatment • High risk (Braden ≤ 12), or • Multiple Stage I's to Single Stage II on trunk or pelvis	Aggressive Treatment or Prevention • High risk (Braden ≤ 12), or • Multiple Stage II's or existence of any deeper ulcers	Critical Treatment or Prevention • High risk (Braden ≤ 12), or • Multiple, more severe ulcers	
PressureGuard Turn Select® E0277 			Can't/won't reposition frequently and has ≤ 2 ulcer-free weight bearing positions or can't be positioned off of ulcer(s)	Manual positioning prohibited by pain or can't be positioned off of ulcer(s) and requires customized turn pattern or >10 minute cycles	Powered Treatment
PressureGuard Easy Air® E0277 			Excessive skin moisture (maceration) exacerbates negative effects of pressure or Can't tolerate repositioning/rotation	Excessive skin moisture (maceration) exacerbates negative effects of pressure or Can't tolerate repositioning/rotation	
PressureGuard® APM/APM² E0277 			Can't/won't reposition frequently and has ≤ 2 ulcer-free weight bearing positions and breakdown does not respond to lower level products	Can be repositioned manually and can be positioned off of ulcer(s)	
PressureGuard CFT® E0373 		Can't/won't reposition frequently and can be positioned off of ulcer(s)	Repositions frequently or infrequently and has ≥ 2 ulcer-free weight bearing positions and can be positioned off of ulcer(s)	Can be repositioned manually and can be positioned off of ulcer(s)	Non-Powered Treatment
PressureGuard Renew® 		Can't/won't reposition frequently and can be positioned off of ulcer(s)			
Geo-Mattress® E0184 	Physical profile indicates need for protection	Reposition frequently by self or caregiver and can be positioned off of ulcer			Non-Powered Prevention
Geo-Matt® Overlay 	Physical profile indicates need for protection	Reposition frequently by self or caregiver and can be positioned off of ulcer			